

Horizon Discovery X-MAN™ cell lines reveal how PI3K mutation can result in EMT switch and increased invasiveness

Sue Griffin¹, Jeffrey Wallin², Jane Guan², Kyle Edgar², Wei Zhou², Ross Francis², Anthony Torres², Peter Haverty³, Jeffrey Eastham-Anderson⁴, Sabrina Arena⁵, Alberto Bardelli⁵, John Goodall¹, Kyla Grimshaw¹, Klaus Hoefflich², Marcia Belvin², Lori Friedman² and Christopher Torrance¹

Horizon Discovery Ltd., Cambridge, UK¹ | Departments of Cancer Signaling and Translational Oncology², Bioinformatics³, Pathology⁴, Genentech, Inc., South San Francisco, CA, USA
Laboratory of Molecular Genetics, Institute for Cancer Research and Treatment, University of Torino Medical School, Candiolo, Italy⁵ | FIRC Institute of Molecular Oncology, Via Adamello 16, Milan, ZIP 20139, Italy⁶

INTRODUCTION


A major challenge facing the development of PI3K-targeted pharmaceutical agents is the lack of patient-relevant, *in vitro* model systems. Studies into the function of oncogenic mutations are commonly carried out in cell line panels that are mutant or wild-type for the gene of interest. Interpretation of these studies can be confounded by additional genetic differences between each cell line background, calling into question the reliability of the results. An alternative approach is to over-express the protein of interest. However, generating a significantly higher protein level than is normally found in the cell can lead to results that do not reflect the true biology of the disease.

Horizon Discovery’s proprietary rAAV GENESIS™ technology directly addresses this problem through the development of X-MAN™ isogenic cell lines. Isogenic cell lines are genetically identical except for the mutation status of an endogenous gene, thereby permitting evaluation of a single genetic alteration expressed at the endogenous level. The technology permits not only knock-out of existing mutations, but also introduction of single or multiple mutations, deletions, amplifications and translocations. Therefore, isogenic cell line models allow the interrogation of protein function without the problems associated with protein over-expression studies or genetically diverse cell line panels.

In order to investigate the effect that a PI3K mutation might have on tumour development, we employed isogenic cell lines derived from a ‘normal’ background (MCF10A mammary epithelial cells). A PI3Kα (H1047R) mutation (commonly found in breast cancer patients), was knocked-in to produce a pair of MCF10A cell lines; parental cells that express wild-type PI3K or PI3Kα (H1047R) cells that express mutant PI3K¹.


CELL LINE GENERATION

Figure 1 shows the generation of the MCF10A PI3Kα (H1047R) isogenic cell line using Horizon Discovery’s GENESIS™ genome editing technology.




RESULTS


Parental and PI3Kα (H1047R) cells were profiled in 2D proliferation assays using a PI3K inhibitor, GDC-0941. Surprisingly, the parental and PI3Kα (H1047R) cell lines displayed similar sensitivities to the inhibitor in this 2D assay format (Figure 2A). Consistent with the principle that PI3K mutations activate AKT and other downstream signalling molecules, the phosphorylation levels of several downstream pathway components were increased in PI3Kα (H1047R) cells in comparison with parental cells (Figure 2B). As expected treatment with the PI3K inhibitor resulted in reduced phosphorylation of the downstream pathway proteins.




Since a PI3K inhibitor did not selectively inhibit the growth of PI3Kα (H1047R) cells, the consequence of introducing the PI3K mutation was further investigated using a microarray study to compare the gene expression profiles of parental and PI3Kα (H1047R) MCF10A cells (Figure 3). This analysis revealed that the PI3Kα (H1047R) cells have up-regulated expression of mesenchymal genes and down-regulated expression of epithelial genes – a profile that is suggestive of the PI3Kα (H1047R) cells undergoing epithelial to mesenchymal transition (EMT).



Increased invasive potential is commonly associated with EMT and oncogenic PI3K is believed to play a critical role in driving metastasis². Therefore, the phenotype of parental and PI3Kα (H1047R) cells grown in 3D was investigated. In this format, PI3Kα (H1047R) cells formed a distinct networked morphology, suggestive of an invasive phenotype, while parental MCF10A cells expressing wild-type PI3Kα grew as small rounded cells (Figure 4). Treatment with the PI3K inhibitor GDC-0941 reversed this morphology but treatment with an AKT inhibitor (AKT1/2i) did not.



To investigate if the morphology displayed when the cells are grown in 3D is linked to increased invasiveness, invasion assays were performed with parental and PI3Kα (H1047R) cells. The PI3Kα (H1047R) cells invaded approximately 5-fold faster than the parental cells over 24h (Figure 5). Furthermore, treatment with GDC-0941 reduced this invasive phenotype, suggesting that the increased invasive capacity is due to the PI3Kα (H1047R) mutation. Interestingly, use of the AKT inhibitor increased the invasive phenotype of the PI3Kα (H1047R) cell line, suggesting that inhibition of AKT can actually result in feedback activation of this pathway.



CONCLUSION

- ❖ By using Horizon’s rAAV GENESIS™ technology, it has been possible to demonstrate that introduction of a PI3K mutation results in EMT switch and increased cell invasiveness.
- ❖ This phenotype was revealed by the ability to engineer expression of the mutation at the endogenous level.

REFERENCES

1. Di Nicolantonio *et al.*, PNAS. 2008. 105(52):20864-9.
2. Samuels *et al.*, Cancer Cell. 2005. 7(6):561-73.